## TAKE AWAY PAIN ...

... is drive, philosophy and task all in one for the therapist Simone Märten. For 20 years now, Simone Märten, masseuse, medical swim advisor, physic therapist for sport, treats versatile illnesses. She has now developed a device which amazes the experts, patients and everybody using it. "PERL" gets to the basis of our feet What is "PERL" capable of - we asked the inventor Simone Märten.

## therapist grow to develop such a technical device?

I have always been fascinated with my high interest in technics on any new development. Of course do I love to work with my hands as a therapist. But that can't be all. Technic has by today grown so much that it has many ways of placing them into therapy. So I combined my affinity for technics with my practical knowledge and experience and feeling what patients need to develop this BMS foot treatment device.

### How does the idea of a physio- How does this health and well- "PERL" has many options to be ness instrument work?

Perl is stimulating the muscles in feet and legs bio-mechanically, it is based on a smooth vibration of 20 Hz. The apparatus moves the feet, encourages the blood circulation, trains the connective tissue and vane flaps, it encourages the drainage of the veins in the legs. For this whole positive result one only needs 5 min. per day. "PERL" stimulates the activity of nerves and therefor achieves a positive result in the mental situation of the user. The device is easy to handle, the weight of 45kg not to heavy and it can be connected to any plug.

used, on diabetics, oncological patients etc. - which ones are the most important to you?

I can't just answer it as there are so many individual needs and complaints. If someone needs a treatment before a long-distance flight against thrombosis, another one not being able to walk by himself anymore. For me "PERL" has proven excellent in results and gone through my tests with full marks. So I know that it helps chemo-patients to ease their side effects, pain and numbness in their feet, it makes it more bearable or takes the pain totally.



Relaxation, mobilization, blood circulation and foot reflex zone stimulation

- bio mechanical muscle stimulation (without pressure points)
- · spends reaxation and comfort
- suitable for hotels, wellness centers, sporting clubs, airport Vip lounges senior homes, fairs, conference centers





# AGIL BODY Gerätevertrieb GmbH

AGIL BODY Gerätevertrieb GmbH Rudolf-Breitscheid-str. 61 14532 Kleinmachnow Mobil: +49 177 24 80 704 info@agil-body.com www.agil-body.com



bilization of your feet is the only one existing in the market in your feet which can occur in cases of oncological sickwhich is based on biochemical muscle stimulation.

mal position of the whole body when using the bio mechanic muscle stimulation (based on 20 Hz) results in a proofed complete mobilization and relaxation.

on the vibration highly effectively through your feet into ents and other patients that suffer from blood circulation

Used on a regular basis a muscle stimulation will achieve a positive influence of remodeling the contours of your legs.

The simple 5 min. treatment with PERL shows with it's totally mechanic effect on connective tissue an impro- Based on the facts and the simple usage we are sure that vement of the veins drainage, training of the valve flaps which herewith result in relief of the cardiac muscle.



This newly developed BMS device for treatment and mo-PERL also offers a pain therapy in the case of neuropathy

It is positioned ergonomically to your feet. With the optiafter chemotherapy on cancer patients.

PERL works effectively with the recommended position of The vibrating area, covered in anti bacterial leather passes legs and feet on diabetics, rheumatics, fibromyalgia patidisorder and all people that want to do something good in their life for their feet.

> Chances and optimal results, if regularly used because of this effective method also offers help in the treatment of balance problems and osteoporosis.

> it will very quickly succeed on the market as a necessary offer in various sectors.

> With you as the presenter of this apparatus with the climbing demand of expectations in the field of wellness and therapy it will surely result in an investment amortizable.

> With this device you attract the attention on to your own

Offers for purchase and financing will be available at our

Simone Märten Managing director www.homaer4you.com

## **Further BMS devices for wellness and training**

#### Profit calculation medical foot BMS-Sigma II treatment device "PERL"

raising expectations in therapy and wellness sectors. We guarantee satisfied customers and their loyalty towards the equipment offered.

To visualize the amortization and lucrativeness in the sector of autonomous self employment as well as an addition to an already existing business or rooms we have calculated a possible income.

The calculated price by us amounts to 5,00 € per 5 min. treatment which does not mean that you can vary that on your own account.

#### **Technical data Parameter Valance** swing frequencies measurements 670 x 670 x 56 mm weight 45 kg

## This results in the following calcu-

2 people using the device 5 min. per hour at an access for 10 hours within 5 days per working week you can earn 100,00 € per day, 500,00 € per week, 2.00,00 € per month without further personnel.

Financing and purchase on a rental basis are available through our bank.

**Technical data** 

With this instrument we meet the The BMS-Sigma II with exchangeable heads / applicators was developed specially for the use on face, eyes, neck, back, extremities and for stimulation of metabolism, nutritive substance of body cells by reaching the extra cellular matrix.

raiailletei	valatice
swing fre-	10-50 Hz
quency	5-50 Hz optional
measurements	200 x 150 x 43 mm
	(control gear BMS-
	Sigma II)
	308 x 62 x 50 mm
	(application part
	BMS-R II)
weight	0,5kg (BMS Sigma II)
	0,8 kg (BMS-R II)

Areas of application		
Swing fre-	Areas to be used at	
quencies		
10-20 Hz	warming up / rege-	
	nerating microcircu-	
	lation	
15-18 Hz	lowering of tonus	
18-22 Hz	anti blocking measu-	
	res	
18-26 Hz	stretching	
20-25 Hz	visceral therapy	
22-24 Hz	improvement of cir-	
	culation, warming up	
22 -26 Hz	muscle toning	
25 Hz	intensify tonus	
from 25 Hz up	mobilization of joints	

#### **BMS-Theta**

#### **Special characteristics**

- maximal comfort and optimal positioning for therapist and patient by ergonomic forms
- · easily reachable operating elements
- very easy control facility
- · easy access to the whole vibrathode as there is no limiting edge
- · Vibrathoden-height and -inclination continuously electronically adjustable
- · vibrathode is suitable for bigger muscle areas as well as for indirect smaller and deeper positioned muscles

Technical data		
Parameter	Valance	
swing fre-	15 - 40 Hz	
quencies	5-50 Hz optional	
measurements	1250 x 570 x 650	
	(1000)	

#### (1090) mm weight 95 kg

Areas of application		
Swing fre-	Areas to be used at	
quencies		
10-20 Hz	warming up/regene-	
	ration/microcircula-	
	tion	
15-18 Hz	lowering of tonus	
18-22 Hz	anti blocking measu-	
	res	
18-26 Hz	stretching	
20-25 Hz	visceral therapy	
22-26 Hz	improvement of cir-	
	culation, warming up	
22-26 Hz	muscle toning	
24-26 Hz	treatment of spasm	
25 Hz	intensify tonus	
from 25 Hz	up mobilization of joints	
27-30 Hz	pain therapy	

On request training/instruction available on BMS devices for therapists and patients.





